**ATEM FLAT FOLD  
DISPOSABLE RESPIRATOR  
EN 149 : 2001  
+A1: 2009**

**INSTRUCTIONS FOR USE :**  
Single Use Filtering face-piece respirators

**IMPORTANT :**  
Single Use Respirators are classified according to one of the categories – FFP1 NR, FFP2 NR, FFP3 NR. Before reading the following information, check the category to which the respirator belongs this is indicated on the packaging and on the respirator.  
**If the mask is marked D, it has passed the dolomite clogging test and suitable for use in very dusty environment. It is the user's responsibility to ensure that the respirator provides the necessary level of protection for the type and concentration of the contaminants (s) in the area where the respirator is intended for use.**

**DIRECTIONS FOR USE :**  
- Fit the respirator and check the leak tightness BEFORE entering the contaminated area  
- Wear the respirator for the whole duration of exposure to contaminants  
- Use the respirator in accordance with applicable health and safety regulation  
- Discard the respirator and replace with the new one if the respirator start to cause breathing difficulty or discomfort or damaged or the smell of vapour present becomes detectable  
- Leave the contaminated area if dizziness, irritation or other distress occurs  
- Only for single use. No maintenance necessary. Do not store and re-use after single use (one shift, maximum usage)  
- Keep un-used respirator in their closed box and store in a dry non-contaminated area between -20 and +40 deg. C. at relative humidity under 75%  
- This mask will only provide protection against the risks and hazards identified in this instruction manual and is suitable for use against solid, water and liquid based aerosols.

**NOTIFIED BODY : BSI Group The Netherlands B.V ( 2797) Say Building John M. Keynesplein 9, 1066 EP Amsterdam**

**CHECKS PRIOR TO USE:**  
Inspect respirator, including any component parts, for holes, tears or damage prior to use. Should any damage be identified, DO NOT USE IT. Get a new one.

**FITTING INSTRUCTIONS :**

1. HOLD RESPIRATOR HORIZONTALLY WITH NOSEPIECE ON TOP AND PULL OPEN ( AS SHOWN)
2. PLACE CHIN AND FACE IN MASK WITH NOSEPIECE ON THE BRIDGE OF YOUR NOSE.
3. PLACE THE TOP BAND HIGH ON THE BACK OF YOUR HEAD AND THE BOTTOM BELOW YOUR EAR
4. ADJUST THE STRAP TENSION BY PULLING THE STRAPS. (FOR THE MASK WITH ADJUSTABLE CLIP ONLY)
5. USING BOTH HANDS ,MOLD METAL NOSEPIECE TO SHAPE OF NOSE FOR TIGHT FIT
6. CHECK FOR LEAK-TIGHTNESS (FACIAL FIT) AS FOLLOWS:
  - a) Place both hands over the respirator and exhale sharply \*positive pressure should be inside the respirator\* if you feel air escaping around the edges ,re-adjust the respirator by tightening the nose-piece and/or increasing tension of the head bands\*repeat the procedure until satisfactory face-fit has been achieved.
  - b) Test the fit for respirator with exhalation valve .Cup both hands over the respirator and inhale vigorously . If no negative pressure felt inside the respirator , re adjust the position of the respirator according to steps #2 to #5 till get a better fit

**IMPORTANT :** In the case of respirators with non-adjustable head bands , if you feel that the head bands are too loose, you can increase tension by tying a small knot in the head band

Please be aware of the appearance of exhalation valve has been changed on ATEM 3132, ATEM 3043 respirator.  
But it won't affect the filtering performance and the valve leakage rate.

**LIMITATION TO USE :** Do not use these respirator or enter or stay in an area where :

- ➔ Oxygen concentration is less than 17% ( UK 19%)
- ➔ Contaminants or their concentrations are unknown or immediately dangerous to life or health
- ➔ Particulate concentration exceed levels fixed by applicable health and safety regulations or protection factor XNPF whichever is lower / Nominal protection factors : 4.5 for FFP1 NR , 12.5 for FFP2 NR and 50 for FFP3NR
- ➔ Gases and/or vapors are present -except for respirators designated as protecting against organic vapours in which cases vapour concentration should not exceed NPF .Respirators protecting against organic vapours should not be used as protection against acid vapours and vice-versa.

**WARNING :**


Do not use for fire -fighting .These respirator do not supply oxygen .Do not use in oxygen deficient atmospheres -eg: tanks or other poor ventilated areas. Do not use in explosive atmospheres. Requirements for leak tightness are unlikely to be achieved if the respirator is worn against beard or facial stubble. When transporting the product , retain in the original packaging and keep away from mechanical and chemical hazards

**IMPORTANT :** If you CANNOT achieve a proper fit, DO NOT enter the contaminated area. See your supervisor for proper fitting respiratory protection

**MANUFACTURER IDENTIFICATION:**

Type: FFP1 NR D ; FFP2 NR D ; FFP3 NR D

NR is to be worn for a single shift only  
D: Meet the Dolomite Clogging test  
CE Standard considered: EN149:2001+A1:2009  
CE Marking : CE 2797

Meaning of pictograms :-  
a) Read instructions before use  
b) Storage conditions  
c) Relative humidity  
d) Expiry date

**Technical information**

| Tests Type                                   | FFP1 NR D  | FFP2 NR D  | FFP3 NR D |
|--|------------|------------|-----------|
| Total Inward Leakage                         | < 22%      | < 8%       | < 2%      |
| Penetration                                  | < 20%      | < 6%       | < 1%      |
| Carbon dioxide content of inhalation air (%) | < 1%       | < 1%       | < 1%      |
| Inhalation resistance<br>30 l/min (mbar)     | < 0.6 mbar | < 0.7 mbar | < 1 mbar  |
| Inhalation resistance<br>95 l/min (mbar)     | < 2.1 mbar | < 2.4 mbar | < 3 mbar  |
| Exhalation resistance<br>160 l/min (mbar)    | < 3 mbar   | < 3 mbar   | < 3 mbar  |

PPE's controlled by :

BSI Group The Netherlands B.V Say Building John M. Keynesplein 9, 1066 EP Amsterdam

Our filtering half masks are in conformity with PPE Regulation (Eu) 2016/425

EU Declaration of Conformity is available to download from our website under Europe of CERTIFICATE

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**ATEM Company Pty Ltd**